

# Electron Phase Microscopy Using Field-Emission Electron Beam

Akira Tonomura<sup>1,2,3</sup>

<sup>1</sup>Advanced Science Institute, RIKEN, Wako, Saitama 351-0198, Japan

<sup>2</sup>Advanced Research Laboratory, Hitachi, Ltd., Hatoyama, Saitama 350-0395, Japan

<sup>3</sup>Funding program for world-leading Innovative Research and development on Science and Technology, Hatoyama, Saitama 350-0395, Japan

Electron phase microscopy has been used to illuminate fundamental phenomena in magnetism and superconductivity by visualizing quantitative magnetic lines of force.

We can use the phase information of field-emission electrons to directly image the quantum world by displaying the phase contours of the transmitted electrons [1]. Projected magnetic lines of force can be quantitatively observed by using holographic interference microscopy [2] based on the Aharonov-Bohm (AB) effect principle [3], and quantized vortices in superconductors can be dynamically observed by Lorentz microscopy. [4]

We experimentally observed the magnetic fields of magnetic heads and superconductors using 1-MV and 300-kV field-emission transmission electron microscopes. An example of tiny magnetic heads is shown in FIG. 1. The head is used for high-density perpendicular recording. It must produce a very high and yet localized magnetic field. However, this is not sufficient. For example, we also need a zero field outside the head in some cases.

In the case of an FeNi head, straight magnetic lines inside the head (FIG. 1 (a)) were changed to “U-turn” lines [5] (FIG. 1 (b)) when the applied magnetic field was reduced from 80 Oe to zero. At the pole tip, the magnetic lines formed a tiny closed circuit, creating vanishing leakage fluxes outside the head. This remanent state is reversible but slightly dependent on the orientation of the previously applied magnetic field (see FIG. 1 (c)). Hysteresis became noticeable when the head material was changed from FeNi to CoNiFe (see FIG. 1 (d)).

We also observed the nucleation and growth of ferromagnetic domains in colossal magnetoresistance materials, in which an electric current begins to flow when the magnetic field increases or the temperature decreases. The microscopic mechanism for this is considered to be due to the fact that isolated tiny ferromagnets nucleate and coalesce in the “insulating” phase, thereby creating conduction paths of electrons. We investigated this directly using Lorentz microscopy and interference microscopy. [6]

The growth of ferromagnetic domains was dynamically observed by Lorentz microscopy (FIG. 2 (a)) while reducing the sample temperature. Interference microscopy revealed that “sepia” leakage flux assisted the domain growth, as shown in FIG. 2 (b).

## References

- [1] A. Tonomura, *Proc. Natl. Acad. Sci. USA*, 102 (2005) 1452.
- [2] A. Tonomura et al., *Phys. Rev. Lett.* 44 (1980) 1430.
- [3] Y. Aharonov and D. Bohm, *Phys. Rev.* 115 (1959) 485.
- [4] K. Harada et al., *Nature* 360 (1992) 51.
- [5] J. J. Kim et al., *Appl. Phys. Lett.* 92 (2008) 162501-1.
- [6] Y. Murakami et al., *Nature Nanotech.* 5 (January 2010) 37.




FIG. 1. Interference micrographs showing magnetic lines of force of magnetic heads for perpendicular recording: (a)  $H = 80$  Oe (NiFe); (b)  $H = 0$  (NiFe) after the upward magnetic field; (c)  $H = 0$  (NiFe) after the downward magnetic field; (d)  $H = 0$  (CoNiFe).




FIG. 2. Growth of ferromagnetic domains in colossal magnetoresistance manganite,  $\text{La}_{0.25}\text{Pr}_{0.375}\text{Ca}_{0.375}\text{MnO}_3$ , thin film at  $T = 62.4$  K: (a) Lorentz micrograph; (b) interference micrograph.